

An Endoscope With 2 DOFs Steering of Coaxial Nd:YAG Laser Beam for Fetal Surgery [Yamanaka *et al.* 2010, IEEE trans on Mechatronics] Automatic Tracking Algorithm in Coaxial Near-infrared Laser Ablation Endoscope for Fetus Surgery [Hu et al. 2014, International Journal of Optomechatronics]

ENGINEERING

GEORGE DWYER

SURGICAL VISION GROUP, CENTRE FOR MEDICAL IMAGE COMPUTING





Clinical Motivation

- Ŵ

- Twin to Twin Transfusion Syndrome
 - Uneven blood flow between twins sharing a placenta
 - High Mortality rate
- Treatment: Fetoscopic Laser Photocoagulation
 - Laser fibre introduced endoscopically
 - Connecting blood vessels in the placenta are ablated

















SURGICAL ROBOT VISION





Current Limitations (Justifications)

- Fibre endoscopes require high power illumination
- Procedure requires high dexterity to target vessels
 - 2D images unknown distance
 - Placenta is delicate

Aims

- Rigid Endoscope to reduce illumination required
- Target the laser without moving the endoscope



GIF









System Framework



SURGICAL ROBOT VISION

5

ROYAL ACADEMY OF ENGINEERING

EPSRC

Engineering and Physical Sciences

wellcometrust







GIF

urq









- ▶ FoV: 70°
- Laser Spot Properties
 - ▷ d = 10 mm, Ø = 2.2 mm
 - > d = 20 mm, ∅ = 3.2 mm
- No illumination channel

450g



7

mic









GIF burg Calibration



Intrinsic Extr

Extrinsic









Results: Bench Testing



Power

- Transmission efficiency = 34%
- Maximum laser power = 14.86 W (50 W)

Area

- 0.9 P_{max}: d = 20 mm; x = 10.7 mm; y = 9.7 mm
- ▶ Firing angles (optical axis)
 - > x: -13 to 11; y: -16 to 6
- Max intensity was out of principal point of mirrors

Accuracy

research group

- Water test; sheet pitch 5mm; target intersections
- Positioning Error(20 mm) = 0.3mm ± 0.2mm
- Error > 1 mm near the edge of the FoV



Power at 20 mm against optical angle

wellco



ROYAL ACADEMY OF ENGINEERING



Results: in vitro Irradiation



cmic

Centre for Medical Image Computi











System Discussion

- Clinically used Systems
 - Rigid Curved Instruments
 - Manually manipulated
 - Diameter < 3mm</p>
 - Built in light source
 - > 2 instrument channels

- Presented System
 - Firing angles of over 20° horizontally and vertically
 - Misaligned
 - Large beam diameter (3.2 mm at 20 mm)
 - Diameter = 7 mm
 - No light source
 - 34% transmission efficiency













System Discussion

Ŵ

- Clinically used Systems
 - **Rigid Curved Instruments**
 - Manually manipulated
 - Diameter < 3mm
 - Built in light source
 - 2 instrument channels

Future System

- Refine manufacturing and placement of optics
- Use of smaller optics to reduce size
- Increase transmission efficiency

12





ROYAL ACADEMY OF







Automatic Tracking Algorithm in Coaxial Near-Infrared Laser Ablation Endoscope for Fetus Surgery











Clinical Motivation and Aim



- Target point movement from respiratory and cardiac motion reduces accuracy
 - Respiratory arrest from general anaesthetic can be applied
 - Postoperative complications
 - For MRI scans, mother hold hers breath (20s)
 - Coagulation takes longer than 20s in total
- Solution

esearch group

Object point tracking algorithm



14

ENGINEERING





Speed is the main consideration



wellcometrust







SURGICAL ROBOT VISION

-=== research group

Preprocessing



- Camera Calibration (Zhang)
 - Remove distortions
 - > Determine relationship
- Reduce image size
- > Apply LoG filter to image
 - > Find edges
 - > Smooth
 - Grey scale

LoG =
$$\frac{x^2 + y^2 - 2\sigma^2}{2\pi\sigma^6} \exp\left(-\frac{x^2 + y^2}{2\sigma^2}\right)$$



EPSRC

Ingineering and Physical Scien

wellcometrust

16

ROYAL ACADEMY OF ENGINEERING



Extracting Keypoints



- Feature from Accelerated Segment Test (FAST)
 - > High computational efficiency
 - Corner detecter
 - Intensity difference to threshold
 - For 768 x 288 image
 - FAST takes 1.5 ms
 - > 50 times lower than SIFT



17







ROYAL ACADEMY OF ENGINEERING



Optical Flow

- Track features from FAST
- Pyramidal Implementation of Lucas-Kanade Optical flow
 - Finds the brightness gradient on a image neighbourhood
 - Mistaken points are removed by Random Sample Consensus Test using a homography matrix (H)



Ŵ













- Homography Matrix
 - Applies relationship between projection of H and relative positions of keypoints in the pre-frame

$$\lambda \begin{pmatrix} q_{cur} \\ 1 \end{pmatrix} = H \begin{pmatrix} q_{pre} \\ 1 \end{pmatrix}$$

- q is object position
- > λ is the z axis when projected in 3D



ROYAL ACADEMY OF

ENGINEERING

Engineering and Physical Science







- Relative position of surrounding points
 - Current position can be calculated from:
 - The current keypoint position s'i
 - Relative position of the keypoint and object position of the preframe

$$q_t = \frac{\sum_i w_g(t_i) w_v(s'_i) (q_{pre} + t_i)}{\sum_i w_g(t_i) w_v(s'_i)}$$

 \triangleright Wg is weight from relative position, Wv is weight from distance to centre of image



research aroup



ACADEMY OF ENGINEERING





- Object position calculated using Homography
- If Homography cannot be counted on switch to relative positions
- Cross correlation coefficient calculated to ensure correct tracking

$$\mathbf{c} = \frac{\sum_{r_t < |q-p| \le s} L_{cur}(p) L_{pre}(p)}{\sqrt{\sum_{r_t < |q-p| \le s} L_{cur}(p)^2} \sqrt{\sum_{r_t < |q-p| \le s} L_{pre}(p)^2}}$$

Kernel size = 2s + 1; L is brightness of LoG image; p is point in kernel

- c > 0.85; object is tracked
- c < 0.85; object is lost; template matching is applied</p>







Experiments

- Respiratory motion compensation
 - Modelled as sinusoidal motion
 - Assumed displacement of 10 mm
 - Maximum Velocity of 10 mm/s
 - > Period of 4s

research group

- Simulate movement by fixing object and moving endoscope
 - Tip displacement measured by a laser displacement sensor





ROYAL ACADEMY OF ENGINEERING





cmic

Centre for Medical Image Comput

		Processing time [ms]	Similarity	Extracted points	Correctly tracked points	Tracking performance
Phantom	10 mm	30.89 ± 2.34 31.51 ± 2.05	0.93 ± 0.05 0.95 ± 0.02	86 ± 23 85 ± 19	81 ± 23 79 + 20	100%
	20 mm	31.31 ± 2.03 31.75 ± 1.75	0.95 ± 0.02 0.96 ± 0.02	85 ± 19 86 ± 24	79 ± 20 79 ± 20	100%
Pig vessel	30 degree 20 mm	31.54 ± 2.57 44.66 ± 12.19	$\begin{array}{c} 0.92 \pm 0.04 \\ 0.91 \pm 0.01 \end{array}$	$\begin{array}{c} 85\pm21\\ 176\pm50 \end{array}$	$\begin{array}{c} 80\pm20\\ 156\pm48\end{array}$	100% 50%

Camera frame rate is 30 fps















wellcome^{trust}







Experiments: under water



- Re-calibrate
- Distance set to 15 mm from centre of object
- Particles
 - > Amniotic fluid particles, 1-2 mm in size at 35 weeks
 - > TTTS is treated before 26 weeks
 - > White paper cut to 1-2 mm is used to simulate particles
- Blurring
 - Blur caused by sudden motion viewed as shift-invariant blurring
 - Simulated using shift-invariant box filter (OPENCV)
 - Blur applied with kernel size of 5 and 11



ENGINEERING











cmic

Centre for Medical Image Computi







ROYAL ACADEMY OF ENGINEERING





		Processing time [ms]	Similarity	Extracted points	Correctly tracked points	Tracking performance
Floating particles		26.31 ± 3.29	0.95 ± 0.01	82 ± 21	73 ± 21	91.39%
Blurring	Original	24.21 ± 2.54	0.96 ± 0.02	85 ± 18	78 ± 21	100%
	Kernel = 5	22.84 ± 1.79	$0.97 \pm .01$	84 ± 23	77 ± 20	99.16%
	Kernel = 11	23.55 ± 2.94	0.97 ± 0.01	97 ± 22	88 ± 24	97.63%



cmic

Centre for Medical Image Computi

ROYAL ACADEMY OF ENGINEERING

EPSRC Engineering and Physical Sciences Research Council

wellcometrust



urg Results: under water





Tracking with Floating Particles



GIF









Blurring





EPSRC Engineering and Physical Sciences Research Council

wellcometrust

Original Image



29

cmic

Centre for Medical Image Computi

ROYAL ACADEMY OF ENGINEERING





cmic

Centre for Medical Image Computi

ROYAL ACADEMY OF ENGINEERING

Blurring





EPSRC Engineering and Physical Sciences Research Council

wellcometrust









Blurring





EPSRC Engineering and Physical Sciences Research Council

wellcometrust

Kernel = 11



cmic

Centre for Medical Image Computi

ROYAL ACADEMY OF ENGINEERING



- FAST & LK Optical Flow
 - Detect and track feature points
 - RANSAC used to remove bad tracking
- Average errors were approximately 10 pixels
- Similarity was on average more than 0.9
- Processing Speed averaged 30 ms
- Future Work Adjust algorithm for better handling of
 - Floating particles
 - Unobvious difference between vessel and background (pig vessel)











